

Explore Tracked Images

Explore Tracked Images allows visualization of the tracked images saved in the output folder as specified by the user. Note: If a **Raw Images Directory** is supplied within the **Folders Tab**, the Cell Tracker will allow contour images to be displayed. Otherwise, only the tracked image masks will be displayed.

The **Show Cell Number(s)** input parameter allows all of the data or only specific cell data to be plotted. The syntax for specifying specific cell numbers to be displayed is a comma separated list. For example, "1,3,5,9" will only display cells "1,3,5,9." A colon can also be used to specify a range of cells. For example, "1,4:8,12" would display cells "1,4,5,6,7,8,12."

The save images button will open a window to a specified directory and save the images as displayed in that directory with a common name as specified in the **common_image_name_** edit box.

To jump to a specific frame, enter the frame number at the end of the slider bar and choose "Go" or hit the "Enter" button on your keyboard. You can also move the slider by clicking with the left mouse button or by using the left and right arrow after clicking at least once on the slider bar with the mouse.